


# Make every drop of water count

Large surface of water, results in quick recovery of humidity in our compact chamber.



## Single Drawer Incubator SD-830

# mini-Cube



Open the "Drawer" to new possibilities

# A clean Unit is a Safe unit



Autoclave-able drawer can be removed entirely. The solid base drawer is made from grade A aluminum and is easy to clean with distinctly designed curved corners.

### SPECIFICATIONS

#### MODEL SD-830

External dimensions (mm)	W310 x D315 x H190
Number of chambers	1
Usable dim. of chamber (mm)	W92.5 x D157 x H25
Operational environment	20~30C
Temp. accuracy	±0.1C
Temp. uniformity	±0.15C
Humidity	95% ±4%RH
Humidity system	Natural evaporation
Required gas	CO <sub>2</sub> /N <sub>2</sub>
Power	AC110-115V 0.8A (Max) 50/60Hz
	AC220-240V 0.6A (Max) 50/60Hz
Weight	15kg
Two Divided door	190 x 127 (×2doors)

